




Product #6638

TechNiche Phase Change Cooling **BACK WRAP**

Powered by CoolPax™

- Secure comfortably and firmly around your back
- Provides cooling relief to aches and pains
- Cool Pax™ Phase Change Cooling inserts maintain a comfortable 58°F (14°C)
- 100% Cotton Carrier
- Includes one back wrap and four (4) 6668 CoolPax™ inserts
- Size: S/M, L/XL
- Color: Blue



POWERED BY

COOLPAX™

CoolPax™ is our Phase Change Material (PCM) cooling technology. It uses a non-toxic liquid that freezes at 58°F/14°C and remains at that temperature for up to 1 hour.

HOW TO USE :

1. Place CoolPax™ inserts in freezer or ice water for 35-40 minutes
2. Insert activated CoolPax™ into garment pockets
3. CoolPax™ vest or accessory is ready to wear!

COOLPAX™ Inserts

Polyurethane Outer Material
Phase Change Material: 58°F
Activation: Ice Water or Refrigerator
Activated Duration: up to 1 hour